



XTENCHO™ RACE TECHNICAL GUIDE

Mammalian Protein Expression
Guide and Protocol

Catalog number PX-XTE-R-001



✉ support-xtencho@proteogenix.fr

☎ +33 (0)3 90 20 54 70

🌐 www.proteogenix.science


ProteoGenix
Accelerate your therapeutic research !



Limited Use Label Licenses

Limited Use Label License for internal Research and Development Use of Mammalian Cell Lines

Notice to Purchaser: The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product and progeny, to perform internal research and development for the sole benefit of the purchaser. The purchase of this product does not grant the purchaser any additional rights, including, without limitation: (a) the right to use the product or its components in commercial applications of any kind, including bioproduction (e.g. use of the product to manufacture therapeutic agents or diagnostic test components, such as proteins, antibodies, viral particles or virus-like particles), quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration; (b) the right to use the product or its components as a therapeutic agent or diagnostic test component; (c) the right to use the product or its components to produce material for use in human clinical trials; or (d) the right to transfer or resell the product or its components in any form, progeny or derivative. For information on obtaining additional rights, please contact us at contact@proteogenix.fr, Licensing and Commercial Supply, ProteoGenix, 19 rue de La Haye 67300 Schiltigheim, France.

TRADEMARKS: All trademarks are the property of ProteoGenix and its subsidiaries unless otherwise specified



Table of Contents

I. Introduction

▶ XtenCHO™ Race Cell Line	03
▶ XtenCHO™ Expression Medium	03
▶ XtenCHO™ Race Transfection Reagent	04
▶ Expression Vector	04
▶ pXten1 Vector Information	05

II. Product Components 06

III. Product Details 07

IV. XtenCHO™ Race Protocol 08

▶ Preparation of working solutions and complete media	08
A. Preparation of complete XtenCHO™ Expression Medium	08
B. Preparation of XtenFect Reagent, Working solution	08
▶ Culture of XtenCHO™ Race Cells	09
A. Recommendations for XtenCHO™ Race Cell culture	09
B. Thaw and establish XtenCHO™ Race Cells	10
C. Subculture XtenCHO™ Race Cells	12
▶ Transfection of XtenCHO™ Race Cells	14
A. Recommendations for transfection of XtenCHO™ Race Cells	14
B. Transfection of XtenCHO™ Race Cells	16

V. Troubleshooting Guide 20

VI. Frequently Asked Questions 23

Ordering Information 32

WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.
NOT FOR HUMAN OR VETERINARY DIAGNOSTICS OR THERAPY.



Part I

Introduction

The XtenCHO™ Race Starter Kit (Cat. no. PX-XTE-R-001) is a complete solution for high-yield transient expression based on suspension-adapted XtenCHO™ Race Cells in a defined, serum-free medium. The XtenCHO™ Race Starter Kit provides cells, expression medium and reagents to transfect 1 liter of cell culture.

■ The XtenCHO™ Race Cell Line

The XtenCHO™ Race Cell line is derived from the Chinese Hamster Ovary CHO-K1 cell line, which has been genetically modified. Combined with the optimized expression vector pXten1 and transfected with the optimized protocol and reagents, XtenCHO™ Race Cells sustain extended transient gene expression, by enhancing plasmid-driven expression and allowing better plasmid maintenance.

The XtenCHO™ Race Cells are adapted to high-density suspension culture in XtenCHO™ Expression Medium. Frozen cells are supplied, and may be thawed directly into XtenCHO™ Expression Medium.

The XtenCHO™ Race Cell line exhibits the following characteristics:

- Derived from the CHO-K1 cell line
- Adapted to high density, serum-free, suspension growth in XtenCHO™ Expression Medium
- Doubling time of approximately 21 hours
- High protein expression

The XtenCHO™ Race Cells are also available separately; see [Page 32](#) for ordering information.

■ XtenCHO™ Expression Medium

The XtenCHO™ Expression Medium is serum-free, animal-component free, and chemically-defined. The XtenCHO™ Expression Medium optimizes productivity with sustained balanced growth and viability post-transfection. The XtenCHO™ Expression Medium supports small- and large-scale transient transfection in XtenCHO™ Race Cells enabling researchers to use the same medium for regular subculturing and for transient transfection.



XtenCHO™ Expression Medium requires additional supplementation of L-Glutamine and Anticlumping agent.

The XtenCHO™ Expression Medium is also available separately; see [Page 32](#) for ordering information.

■ XtenCHO™ Race Transfection Reagent

The XtenFect Reagent is optimized for the transfection of nucleic acids into high-density XtenCHO™ Race Cell cultures grown in XtenCHO™ Expression Medium. The XtenFect Reagent is provided as a 10 X concentrated stock solution, and the XtenFect Reagent Working solution has to be reconstituted prior to use on XtenCHO™ Race Cells. Then the XtenFect Reagent Working solution can be added directly to cells in XtenCHO™ Expression Medium. It is not necessary to pre-form complexes with plasmid DNA before addition to the cell cultures.

***Note:** Use of transfection reagents other than the XtenFect Reagent to transfect high density XtenCHO™ Race cultures can lead to substantially reduced performance.*


The XtenCHO™ Race Enhancer is designed to maintain high viability in high density transfected XtenCHO™ Race Cells, resulting in extended transient gene expression and enhancement of protein production.

■ Expression Vector

The pXten1 Expression Vector is an empty vector optimized for the expression of recombinant proteins or antibodies in XtenCHO™ Race Cells. The plasmid is designed to be used in XtenCHO™ Race Cells permitting enhanced plasmid-driven expression and better plasmid maintenance, resulting in increased levels of protein expression.


The Antibody Expression Positive Control Vector is provided as a control for transfection and expression in XtenCHO™ Race Cells. The control consists of pXten1 plasmids expressing the heavy and light chains of a human IgG1 antibody. The control is a transfection-grade plasmid provided at a concentration of 1 µg/µL with a 1 : 1 light chain / heavy chain ratio and is enough to transfect up to 150 mL of XtenCHO™ Race Cells.

■ pXten1 Vector Information



Sequencing **primers** (after subcloning into the Multiple Cloning Site, **MCS**):

Forward primer (5' - 3')	TCCACAGGTGTCCACTC
Reverse primer (5' - 3')	AACAGATGGCTGGCAACTA





Part II

Product Components

Each XtenCHO™ Race Starter Kit includes the following components:

Component	Quantity/Size	Storage
XtenCHO™ Race Cells (10 ⁷ cells/mL)	2 x 1 mL	Liquid nitrogen ¹
XtenCHO™ Expression Medium	1 L	4°C, protect from light
XtenCHO™ Race Transfection Kit 1 L, contains:	1 kit	
XtenFect Reagent, Stock Solution (dilute prior to use)	1 x 0.51 mL	-80°C, until dilution
XtenFect Reagent, Working solution (reconstitute prior to use)	3 x 1.44 mL	-20°C, until reconstitution
XtenCHO™ Race Enhancer	3 x 1.2 mL	-20°C
pXten1 Expression Vector (1 µg/µL in T10E1 buffer, pH8)	1 x 10 µL	-20°C
Antibody Expression Positive Control Vector (1 µg/µL in T10E1 buffer, pH8)	1 x 150 µL	-20°C

Reagents to be supplied by the user:

- **L-Glutamine** solution, culture grade, 200 mM
- **Anticlumping agent**, culture grade
- **Endotoxin-free**, transfection-grade **plasmids** expressing your protein of interest
- **DMSO**



Part III

Product Details

Kit Size

The XtenCHO™ Race Starter Kit provides cells, culture medium, and reagents to transfect **1 liter** of cell culture.

Shipping Conditions

Dry ice for cells and reagents and blue ice for the medium.

Storage Conditions

Stable for **12 months** from the date of manufacture, if cells, reagents and medium are conserved at their optimal storage temperature. Diluted XtenFect Reagent, Working solution, and XtenCHO™ Race Enhancer should be used within **4 weeks** after dilution and thawing, respectively.

Validation data

Tested with Antibody Expression Positive Control Vector.

Expected Control Results

Productivity > **250 mg/L** (before purification).



Part IV

XtenCHO™ Race Protocol

Preparation of Working Solutions and Complete Media

A. Preparation of complete XtenCHO™ Expression Medium

The XtenCHO™ Expression Medium is formulated without **L-Glutamine** and **Anticlumping agent**. L-Glutamine is quickly hydrolyzed in the medium, decreasing its actual concentration and generating by-products of hydrolysis which can be toxic to the cells.

To ensure optimal L-Glutamine concentration in medium, end users should supplement the XtenCHO™ Expression Medium with the recommended concentration of fresh L-Glutamine (8 mM). An Anticlumping agent (such as 0010057AE from Gibco) should also be added to the cultures at a final concentration of 0.5% for cell subculturing prior to transfection and post transfection.

However, the Anticlumping agent should not be added to the cultures during the transfection or it will result in decreased protein yields. It is highly recommended to remove the Anticlumping agent during the last passage before transfection, by replacing completely the culture medium with XtenCHO™ Expression Medium supplemented with 8 mM L-Glutamine and without Anticlumping agent.

To prepare complete XtenCHO™ Expression Medium:

- Supplement XtenCHO™ Expression Medium with **8 mM L-Glutamine**
- Add **0.5% Anticlumping agent** to the culture when required (cell subcultures and maintenance, and cultures post-transfection)

B. Preparation of Xtenfect Reagent, Working solution

To ensure optimal and prolonged performance of the transfection reagent, the XtenFect Reagent is supplied as a 10 X concentrated stock solution and the XtenFect Reagent Working solution has to be reconstituted prior to use on XtenCHO™ Race Cells. XtenFect Reagent Stock solution should not be used non-diluted on cells or it will dramatically affect cell viability.



To reconstitute the XtenFect Reagent Working solution:

1. Thaw one tube labeled “XtenFect Reagent, Working solution” and let it reach room temperature. The tube contains 1440 μL of XtenFect Dilution Buffer.
2. Thaw the tube of labeled “XtenFect Reagent, Stock solution”. The tube contains 510 μL of 10 X concentrated XtenFect Reagent.
3. Pipet 160 μL of “XtenFect Reagent, Stock solution” and add them to the XtenFect Dilution Buffer into the tube labeled “XtenFect Reagent, Working solution”.
4. Mix by inversion.
5. The reconstituted XtenFect Reagent Working solution is now ready to use.
6. Concentrated XtenFect Reagent Stock solution should be frozen again and stored at -80°C . Avoid multiple freeze-thaw cycles as it can lead to reduced transfection performances. If needed, the XtenFect Reagent Stock solution can be aliquoted.
7. The diluted XtenFect Reagent Working solution could be stored at 4°C and should be used within 4 weeks for optimal performance.

Culture of XtenCHO™ Race Cells

A. Recommendations for XtenCHO™ Race Cell culture

Follow the general guidelines below to thaw, subculture and expand XtenCHO™ Race Cells. All cell cultures must be performed in **microbiological safety cabinet** using aseptic technique to ensure sterility. All media and reagents that come in contact with the cells must be sterile.

Store the frozen cells in liquid nitrogen in the vapor phase until thawing. **Do not store the cells at -80°C** . For all cell manipulations, mix the cells by gentle **swirling**; **avoid extensive shaking/pipetting**. Allow freshly thawed cells to recover in culture for **3 passages** before transfecting.

For maintenance of cells, passage XtenCHO™ Race Cells when they reach a density of approximately **$1.5 \times 10^6 - 2.5 \times 10^6$ viable cells/mL** (i.e., early log-phase growth), typically every **2–3 days**. Please refer to **Table 1** on [Page 12](#) to consult recommended cell seeding densities.



Note: Cells that are subcultured at densities outside of this early log-phase growth window may show longer doubling times and lower titers over time. Modify the initial seeding density to reach the target cell density of $1.5 \times 10^6 - 2.5 \times 10^6$ viable cells/mL at the time of subculturing.

XtenCHO™ Race Cells should be cultured between **2 and 20 passages**. Optimal performances are obtained in these conditions. Use a **hemocytometer** with the trypan blue exclusion method or an **automated cell counter** to determine cell viability. Log-phase cultures should be >95% viable. When thawing cells, **transfer cells into pre-warmed medium**.

B. Thaw and establish XtenCHO™ Race Cells

INTRODUCTION

The protocol below describes steps to thaw the XtenCHO™ Race Cells and initiate cell culture. The XtenCHO™ Race Cells are supplied in a vial containing 1 mL of cells at 1×10^7 viable cells/mL in 90% complete XtenCHO™ Expression Medium and 10% DMSO. Thaw the cells directly into **XtenCHO™ Expression Medium, pre-warmed at 37°C**.

REQUIRED MATERIALS

- XtenCHO™ Race Cells, frozen cryovial
- XtenCHO™ Expression Medium, pre-warmed to 37°C
- L-Glutamine, 200 mM
- Anticlumping agent
- 125-mL disposable, sterile, vented, baffled Erlenmeyer shake flask
- Reagents and equipment to determine viable cell density and viability (e.g., hemocytometer or an automated cell counter, trypan blue)
- Incubator shaker set at 37°C with $\geq 80\%$ relative humidity and 5% CO₂

Note: Store the frozen cells in liquid nitrogen in vapor phase, until thawing. Do not store the cells at -80°C . For all cell manipulations, mix the cells by gentle swirling; avoid extensive shaking/pipetting.



THAWING OF XTENCHO™ RACE CELLS

1. Remove the vial of cells from liquid nitrogen and immediately swirl in a **37°C** water bath for **90 seconds maximum** to thaw the cells rapidly until only a small amount of ice remains. Do not submerge the vial in water.
2. Just before the cells are completely thawed, decontaminate the vial by wiping it with **70% ethanol** before opening it in a laminar flow hood.
3. Quickly transfer the entire contents of the cryovial into **8 mL** of pre-warmed XtenCHO™ Expression Medium. In order to recover all the XtenCHO™ Race Cells, we advise rinsing the cryovial with XtenCHO™ Expression Medium. Mix by inversion.
4. Centrifuge the cellular suspension for **5 minutes at 300 x g**. Discard the supernatant. This step is intended to remove the freezing medium.
5. Resuspend the cells in a small volume of XtenCHO™ Expression Medium supplemented with 8 mM L-Glutamine (**2 mL**).
6. Determine viable cell number and viability.
7. Transfer the total amount of cells to pre-warmed XtenCHO™ Expression Medium supplemented with 8 mM L-Glutamine in a baffled shake flask at a final seeding density of **0.3×10^6 cells/mL**. Add Anticlumping agent at a final concentration of **0.5%**.
8. Incubate the flask in a **37°C** incubator with **≥80%** relative humidity and **5% CO₂** on an orbital shaker platform until cultures reach a density of **$1.5 - 2.5 \times 10^6$ viable cells/mL**.

*Note: Set the shake speed to 140-150 rpm for shakers with a 19 mm shaking diameter, 125-130 rpm for shakers with a 25 mm shaking diameter and 95-100 rpm for shakers with a 50 mm shaking diameter. Please refer to **Table 2** on [Page 13](#) for shaking speed recommendations.*

9. On **Day 2 or 3** post-thaw, determine viable cell density and viability percentage. Cell viability should be **≥90%** by 3 days post-thaw.
10. Continue to monitor cell density and viability and subculture the cells once the culture has reached **$1.5 - 2.5 \times 10^6$ viable cells/mL** (typically 2-3 days post-thaw) using the procedure described below.



C. Subculture XtenCHO™ Race Cells

INTRODUCTION

XtenCHO™ Race Cells are capable of achieving high cell densities; therefore, we recommend that the cells reach a density of **1.5 – 2.5 × 10⁶ viable cells/mL** at the time of subculturing.

Note: Do not allow the cells to reach a density of 3.5 – 4 × 10⁶ and higher during cell subculture.

REQUIRED MATERIALS

- XtenCHO™ Race Cell cultures at 1.5 – 2.5 × 10⁶ viable cells/mL
- XtenCHO™ Expression Medium supplemented with 8 mM L-Glutamine
- Anticlumping agent
- Disposable, sterile, vented, baffled Erlenmeyer shake flask
- Reagents and equipment to determine viable cell density and viability (hemocytometer or an automated cell counter, trypan blue)
- Orbital shaker in a 37°C incubator with ≥80% relative humidity and 5% CO₂
- DMSO

PASSAGE OF XTENCHO™ RACE CELLS

1. Determine the viable cell density of the culture and calculate the volume of cell suspension required to seed a new shake flask according to the recommended seeding densities in **Table 1** and the recommended culture volumes in **Table 2**.

Table 1. Recommended seeding densities for routine subculturing.

Subculture timing	Recommended seeding density
For cells ready 2 days post-subculture	0.3 × 10 ⁶ viable cells/mL
For cells ready 3 days post-subculture	0.2 – 0.3 × 10 ⁶ viable cells/mL



Table 2. Recommended volumes and shaking speeds for routine cell culture maintenance.

Flask size	24 Deepwell plate	125 mL	250 mL	500 mL	1000 mL	2000 mL	3000 mL
Culture volume (mL)	3.5	30–40	60–100	125–200	250–400	500–800	750–1200
Shaking speed							
19–mm orbit	275 ± 5 rpm			140 ± 10 rpm			95 ± 5 rpm
25–mm orbit	240 ± 5 rpm			130 ± 5 rpm			90 ± 5 rpm
50–mm orbit	x			95 ± 5 rpm			85 ± 5 rpm
Vessel type	Deepwell plate	Flask, vented and baffled					

Note: We recommend using 24 Deepwell plates from Axygen (Axygen™ Storage Microplate PDW10ML24CS) as well as sealing films from AeraSeal (ref A9224-50EA).

- Transfer the calculated volume of cells to a new tube and centrifuge for **5 min at 300 x g**. Discard the supernatant. Resuspend the pellet in fresh XtenCHO™ Expression Medium supplemented with **8 mM L-Glutamine** and transfer the cellular suspension to a baffled shake flask. Add Anticlumping agent to the cultures at a final concentration of **0.5%**.
- Incubate flasks in a **37°C** incubator with **≥80%** relative humidity and **5% CO₂** on an orbital shaker platform until cultures reach a density of **1.5 – 2.5 × 10⁶ viable cells/mL**.
- Repeat Steps 1–3 to amplify the cells for transfection.

CRYOPRESERVATION

- Cultivate cells during **3 passages** before freezing, let the cells reach a viable cell density of **1.5 – 2.5 × 10⁶ cells/mL** and **>95% viability** before freezing. Do not freeze cells at <90% viability.
- Determine the viable cell density of the culture and calculate the volume of cell suspension required to freeze aliquots of XtenCHO™ Race Cells, at a density of **1 × 10⁷ viable cells/mL**.
- Calculate the required volume of freezing medium, prepare it and maintain it at 4°C. Freezing medium should be composed of **90 % XtenCHO™ Expression Medium** (supplemented with L-Glutamine and Anticlumping agent) and **10% DMSO**.



4. Centrifuge the cells at **300 × g for 5 minutes**, discard the medium, and replace it with chilled freezing medium. Gently resuspend the cell pellet by pipetting.
5. Aliquot **1 mL** of the cell suspension into each cryovial.
6. Place cryovials quickly in a freezing container designed to achieve a cooling rate of **-1°C/minute**, the optimal rate for cell preservation, then place the freezing container in a -80°C freezer.
7. Transfer frozen vials to liquid nitrogen for long-term storage.

Note: Cells should not be in contact with the freezing medium more than 10 minutes before putting them in the -80°C freezer. Moreover, it is important to achieve a cooling rate of -1°C/minute. When not followed, this may increase cell death during freezing.

Transfection of XtenCHO™ Race Cells

A. Recommendations for transfection of XtenCHO™ Race Cells

INTRODUCTION


Use freshly thawed cells and allow them to recover in culture for **3 or more passages** before transfecting. We recommend not using XtenCHO™ Race Cells aged more than 20 passages for transfection to keep optimal transfection efficiency and yields.

During all cell manipulations, mix the cells by **gentle swirling; avoid vigorous mixing/pipetting**. Cell health is critical to maximal performance.

Make sure you have **replaced** the culture medium by **medium without Anticlumping agent** the day before the transfection.

For optimal transfection of high-density suspension XtenCHO™ Race cultures, use the XtenFect Reagent included in the transfection kit. Make sure to dilute the XtenFect Reagent, Stock solution into the XtenFect Reagent, Working solution before use (**see protocol for preparation of XtenFect Reagent Working solution on [Page 8](#)**). The use of transfection reagents other than the XtenFect Reagent to transfect high density XtenCHO™ Race cultures can lead to substantially reduced performance.

The **plasmid DNA of interest** or **Antibody Expression Positive Control Vector** should be added directly on the cells and cultures, and mixed by gentle swirling.



For expression of recombinant antibodies, we observed that a **1 : 1 mass ratio of heavy and light chain encoding plasmids can be used with the XtenCHO™ Race system**. Most of the cells will be co-transfected with the two plasmids.

The expression of the recombinant antibody will then be dependent on the expression rate of each chain and can be further **optimized**, if necessary, by **modifying heavy to light chain plasmid ratio**. We recommend cloning the heavy and light-chain subunits separately into the **pXten1 expression plasmid**, and then optimizing the ratios of the two plasmids.

It is important to add XtenFect Reagent Working solution **slowly, drop-by-drop** to the cells while swirling the flask. The cell cultures have to be placed in a **37°C** incubator with an orbital shaker, **≥80%** relative humidity, and **5% CO₂**.

The XtenCHO™ Race Enhancer should be added **2 hours post-transfection** and the cell cultures placed back in a **37°C** incubator, under the same conditions.

Unlike some other serum-free media formulations, XtenCHO™ Expression Medium **does not inhibit transfection**. XtenCHO™ Expression Medium is specifically formulated to enable transfection without the need to add feed.

The harvesting time point depends on the nature of the protein expressed and should be adjusted accordingly. Expression levels can vary greatly from protein to protein, and the yield will also be impacted by the stability of the protein.

Transfected cells can be kept until viability drops under **50%** or till **Day 15** post-transfection or **Day 7–8 post-transfection for transfection in 24 Deepwell plates**. Maintaining high cell viability at the time of harvest facilitates downstream processes, such as protein purification.

For stable proteins, as recombinant antibodies, time of harvest occurs generally on **Day 13 to Day 15**. For intracellular or unstable proteins, it might be necessary to harvest the cultures at an earlier time point, to avoid toxicity or degradation. Optimal time of harvest should be determined empirically.

REQUIRED MATERIALS

- XtenCHO™ Race Cells cultured in XtenCHO™ Expression Medium + 8 mM L-Glutamine
- High quality and purity; endotoxin-free plasmid DNA preparations of expression vectors
- Antibody Expression Positive Control Vector
- XtenFect Reagent Working solution, reconstituted from Stock solution, pre-warmed to room temperature (**see protocol for Preparation of XtenFect Reagent Working solution on [Page 8](#)**)



- XtenCHO™ Expression Medium, pre-warmed to room temperature, supplemented with 8 mM L-Glutamine
- Anticlumping agent
- XtenCHO™ Race Enhancer, pre-warmed to room temperature
- Disposable, sterile Erlenmeyer baffled flasks
- Orbital shaker in a 37°C incubator with ≥80% relative humidity and 5% CO₂
- Reagents and equipment to determine viable cell density and percent viability

SCALE-UP OF TRANSFECTIONS

The XtenCHO™ Race Starter Kit allows scaling up and down transfection volumes, from 3.5 mL to 1000 mL. Follow the guidelines indicated in **Table 3** for cell number, transfection and reagent volumes. For large flasks (Fernbach), the shaking speed of the cultures must be adjusted, please see **Table 2** on [Page 13](#) for guidelines.

Table 3. Guidelines for cell number and volumes for transfection at different scales.

Cell culture vessel	Cell number (x10 ⁶)	DNA quantity (µg)	XtenFect Reagent, Working Solution (µL)	XtenCHO™ Race Enhancer (µL)	Transfection volume (mL)	Final transfection volume (mL)
24 Deepwell plate	8.75	5.6	16.8	11.2	1.75	3.5
125 mL flask	75	48	144	96	15	30
250 mL flask	200	128	384	256	40	80
500 mL flask	500	320	960	640	100	200
1000 mL flask	625	400	1200	800	125	250
2000 mL flask	1250	800	2400	1600	250	500
3000 mL flask	2500	1600	4800	3200	500	1000

B. Transfection of XtenCHO™ Race Cells

Determine first the volume of cells you wish to transfect and refer to **Table 3** above, for cell number and corresponding volumes for plasmids and transfection reagents.

1. Subculture and expand XtenCHO™ Race Cells until reaching a density of approximately **1.5 – 2.5 × 10⁶ viable cells/mL**.



DAY -1: PASSAGE OF XTENCHO™ RACE CELLS

2. On the day prior to transfection (Day -1), passage the XtenCHO™ Race culture to a final density of 1.2×10^6 viable cells/mL in XtenCHO™ Expression medium supplemented with 8 mM L-Glutamine. **Do not add the Anticlumping agent to the culture.** Allow the cells to grow for 24 hours.

DAY 0: TRANSFECTION OF XTENCHO™ RACE CELLS

3. On the next day (**Day 0**), determine viable cell density and viability. The cells should have attained a density of approximately $2 - 3 \times 10^6$ viable cells/mL. Viability should be **>90%** to proceed with transfection. Centrifuge, then resuspend XtenCHO™ Race Cells at a final density of 5×10^6 viable cells/mL in half of the final volume wanted with fresh XtenCHO™ Expression Medium supplemented with 8 mM L-Glutamine, following the guidelines indicated in **Table 3** on [Page 16](#). **Do not add the Anticlumping agent**, as it will dramatically impact your transfection efficiency. Swirl the flasks or the plates gently to mix the cells. *For example, for a final 30 mL transfection volume, prepare a cell culture with 5×10^6 cells/mL in 15 mL XtenCHO™ Expression Medium with 8 mM L-Glutamine.*

Note: We recommend discarding the remaining cells and not using high-density cells for routine subculturing.

4. Prepare the XtenFect Reagent Working solution as indicated on [Page 8](#) and allow it to reach room temperature. Mix the XtenFect Reagent Working solution by inversion.
5. Thaw the expression plasmids and mix them by inverting the tubes.
6. Add the expression plasmids directly on the cells in a row. Refer to **Table 3** on [Page 16](#) for DNA quantity. Gently swirl the flasks or the plates to mix. *For example, for a final 30 mL transfection volume, add 48 µg of expression plasmids to the cells.*
7. Add the XtenFect Reagent Working solution slowly, drop-by-drop on the cells, while swirling the flask during addition. Refer to **Table 3** on [Page 16](#) for XtenFect Reagent Working solution volume. *For example, for a final 30 mL transfection volume, add 144 µl of XtenFect Reagent Working solution to the cells.*
8. Incubate the cells in a **37°C** incubator with a humidified atmosphere containing **5% CO₂** in air on an orbital shaker.



DAY 0, 2 HOURS POST-TRANSFECTION: ADD EXPRESSION MEDIUM AND XTENCHO™ RACE ENHANCER

9. Allow XtenCHO™ Expression Medium supplemented with 8 mM L-Glutamine and XtenCHO™ Race Enhancer to reach room temperature. **2 hours post-transfection**, add the other half volume of XtenCHO™ Expression Medium supplemented with 8 mM L-Glutamine to the culture. Add the required volume of XtenCHO™ Race Enhancer directly on the cells. Refer to **Table 3** on [Page 16](#) for XtenCHO™ Race Enhancer volume. Gently mix by swirling the flask. Viable cell density and viability can be checked at this point.

For example, for a final 30 mL transfection volume, add 15 mL of XtenCHO™ Expression Medium (supplemented with 8 mM L-Glutamine) and add 96 µL of XtenCHO™ Race Enhancer to the transfected cells.

Note: *It is normal to observe a drop in cell viability. However, cell viability should not be <80% on Day 0.*

10. Return the flasks or the plates to the **37° C** incubator with a humidified atmosphere containing **5% CO₂** on an orbital shaker.

DAY 1: ADD ANTICLUMPING AGENT AND SHIFT CULTURES TO 33°C

11. **24 hours** post-transfection, add Anticlumping agent to the flask at a final concentration of **0.5%**, gently swirling the flask during the process. Shift the transfected XtenCHO Race Cells to **33°C** incubator with a humidified atmosphere containing **5% CO₂** with shaking. Viable cell density, viability and transfection efficiency can be checked 24 hours and 48 hours after transfection.

Note: *It is normal to observe a drop in cell viability. However, cell viability should not be <60% on Day 1 and Day 2.*

DAYS 10 TO 15: MONITOR CULTURE VIABILITY AND HARVEST WHEN NECESSARY

12. Monitor cell density and viability regularly in the cultures and maintain transfected cells culture until viability drops under **50%** or till **Day 15** post-transfection or **Day 7-8 post-transfection for 24 Deepwell plates**. Note that if protein is not stable or toxic, the time of harvest has to be adapted. When expressing a protein for the first time, you could harvest cells or media at several time points post-transfection to optimize the length of the expression run.



13. Harvest by centrifugation cultures for **5 min at 300 x g** to pellet cells. For recombinant proteins expected to be expressed intracellularly, keep and freeze cell pellets at **-80°C** for further analysis. Culture supernatants can be further clarified by centrifugation for **30 min at 5000 x g** and/or filtration using a **0.22 µm membrane**. It is recommended to proceed quickly with the purification of the protein of interest, or to freeze supernatants at **-80°C** for later use.

TRANSFECTION AND EXPRESSION CONTROL

The Antibody Expression Positive Control Vector is provided as a positive control for transfection and expression in XtenCHO™ Race Cells. The control contains **pXten1 plasmids** expressing the heavy and light chains of a **human IgG1 antibody**. The control is a transfection-grade plasmid provided at a concentration of **1 µg/µL with a 1 : 1 light chain / heavy chain ratio** and is sufficient to transfect up to **150 mL** of XtenCHO™ Race Cells.

The antibody will be expressed and secreted into XtenCHO™ Expression Medium, with optimal yields obtained between **Days 13–15** post-transfection. Typically, yields in crude cell culture supernatants reach concentrations of **> 250 mg/L**, as determined with a Sartorius Octet™ instrument with a protein A biosensor. We recommend to transfect **30 mL of XtenCHO™ Race Cells using 48 µL of the Antibody Expression Positive Control Vector** with the optimal transfection protocol and reagents as described above **or 5.6 µL of the Antibody Expression Positive Control Vector for the 24 Deepwell plates protocol**.



Part V

Troubleshooting Guide

Problem	Solution
---------	----------

Cells show poor growth and viability post thaw. What can I do?

The viability of XtenCHO™ Race Cells can be monitored on the day of thawing and should be >90%. The viability should rise quickly after one or two passages post-thaw and reach >95%.

We advise being very quick when thawing: do not let the cells thaw at room temperature and do not let them stay in the freezing medium for too long time before diluting them with XtenCHO™ Expression Medium.


We advise using freshly supplemented XtenCHO™ Expression Medium: an old preparation is not optimal for cells after thawing as significant degradation of L-Glutamine can occur. Avoid using medium supplemented with L-Glutamine more than 7 days ago. It is also suggested to use pre-warmed complete XtenCHO™ Expression Medium for thawing.

During all cell manipulations, **we recommend avoiding extensive mixing and pipetting of the cells.** It is highly recommended to keep a separate flask for cell culture maintenance and never use cells that have reached high density for routine culture.

For optimal performance, it is critical that the shaking diameter, shaking speed, flask size/type and volume of culture to be transfected **correspond to the recommended conditions** described in the present guide, for both regular subculture and expression experiments.

We advise ensuring the equipment is calibrated for **CO₂**, temperature and humidity. **CO₂ levels should not exceed 6%.**

If poor viability and growth persist, we recommend checking the temperature, CO₂, humidity and shaking of your incubator. Discard cells if you notice morphological changes in the culture. It is highly recommended to thaw another vial of cells rather than trying to keep culturing these cells.



What can I do to avoid clumping of my XtenCHO™ Race Cells?

First, make sure that the XtenCHO™ Expression Medium is supplemented with the **Anticlumping agent**. Then transfer the cells into a new vented and baffled flask if you notice the occurrence of clumps around the neck during subculturing. Clumps may result in the underestimation of the cell density when seeding, leading to overgrown cultures. Such clumps can be also observed after transfection but will not dramatically affect the culture.

Do not try to pick up cells attached around the neck of the flask.

During subculturing, if cell clumps are noticed, resuspend cells carefully and at each passage, allow cell clumps to settle down in the vessel for a few seconds and harvest the suspended cells for subculturing.

Clumping should not affect cell viability, but it can affect growth as cell density cannot be measured accurately if clumps are present in the suspension. We recommend verifying that the cells are cultured according to the recommendations in the present guide in the presence of **0.5% of Anticlumping agent**, below passage 20, and free of contamination.

What can I do if I observe a premature drop in viability post transfection?

Depending on the nature of the protein to be expressed and the expression strategy (**intracellular or secreted expression**), the proper time of harvest should be adjusted. Unstable proteins or toxic proteins can lead to a premature drop in viability.

A drop in viability could also result from suboptimal culture conditions. The shaking speed of the incubator, as well as temperature, CO₂ and humidity levels should be verified. Volume and size of the flask should be appropriate. Depending on the flask brand and the model of the incubator, shaking speed might have to be adjusted.



The expression of my recombinant protein is lower than expected with the XtenCHO™ Race Expression system. What should I do to improve my results?

Combined with the optimized expression vector, pXten1, and transfected with the optimized protocol and reagents, **XtenCHO™ Race Cells** sustain extended transient gene expression, by enhancing plasmid-driven expression and allowing better plasmid maintenance.

Therefore, **we recommend using the pXten1 Expression Vector to express your protein of interest. XtenCHO™ Race Cells** can sustain high expression levels from various expression plasmids but should show optimal performance when using pXten1 as the expression vector.

It is recommended to use **high quality endotoxin-free plasmid** preparations for transfection, as low quality and low purity plasmids will result in low transfection efficiency and increased cell toxicity, dramatically affecting bioproduction yield.

Expression levels can vary greatly from protein to protein; and yield will also be impacted by protein stability. We strongly recommend expressing the positive antibody control to determine if the **low yield** is rather due to low expression of the protein of interest, suboptimal culture conditions or suboptimal transfection conditions. If you are not achieving the expected yield with the antibody positive control, we recommend checking the following:

- Monitor cell viability and verify that XtenCHO™ Race Cells viability is > **95% during subculture** and at time of transfection
- Verify that the doubling time of XtenCHO™ Race Cells is close to **24 hours**
- Verify that your shake speed is about 120–130 rpm for 25 mm orbit shakers and about 140–150 rpm for 19 mm orbit shakers

Further culture optimization may be required.



Part VI

Frequently Asked Questions

1. How does the XtenCHO™ Race Expression System work?

Combined with the optimized pXten1 Expression Vector and transfected with the optimized protocol and reagents, XtenCHO™ Race Cells sustain extended transient gene expression, by enhancing plasmid-driven expression and allowing better plasmid maintenance.

2. How should the XtenCHO™ Race Cells be stored upon receipt?

Upon receipt of the cells on dry ice, place the vials immediately into liquid nitrogen until thawing; **do not store cells at -80°C**. Once thawed and transferred into medium in a vented, baffled shake flask, cells should be cultured at **37°C with 5% CO₂** on a shaker platform set to 125–130 rpm for a shaker with a 25 mm orbital diameter or 140–150 rpm for a 19 mm orbital diameter. Cells should show a high viability upon thawing and should recover quickly post-thaw, reaching their normal 21-hour doubling time within 1–2 passages.

If viability or growth after 1–2 passages are low, cells may have been compromised during either storage, or thawing, or shipment, and should not be used.

3. How should reagents of the XtenCHO™ Race kits be stored upon receipt?

XtenCHO™ Race Cells should be placed immediately **in liquid nitrogen, in the vapor phase**. The XtenCHO™ Expression Medium (1L) should be kept at 4°C. The XtenFect Reagent, Stock solution should be stored at -80°C upon receipt. The **XtenCHO™ Race Enhancer** and the non-reconstituted XtenFect Reagent, Working solution can be stored at -20°C. Plasmids (pXten1 Expression Vector and Antibody Expression Positive Control Vector) can be stored at -20°C.



4. How should reagents of the XtenCHO™ Race kits be handled, prepared and stored?

All reagents need to be brought to room temperature prior to use for transfection. Allow the XtenCHO™ Race Enhancer to thaw and equilibrate at room temperature prior to use. The XtenCHO™ Race Enhancer is provided as a **ready to use solution in 3 aliquots**. The XtenCHO™ Race Enhancer **should not be frozen again**. It can be **stored at 4°C** and should be used within **4 weeks**.

To ensure long term stability of the transfection reagent, the XtenFect Reagent, Stock solution is provided as a **10 X concentrated stock solution**. It must be **diluted 10 times** prior to use and should not be used undiluted on cells.


On the day of transfection, allow the XtenFect Reagent, Stock solution and the non reconstituted XtenFect Reagent, Working solution to thaw and equilibrate at room temperature. Then pipet **160 µL** of XtenFect Reagent, Stock solution into the XtenFect Reagent, Working solution, which contains **1.44 mL** of XtenFect Dilution Buffer, and mix by inversion.

Concentrated XtenFect Reagent, Stock solution **should be frozen again and stored at -80°C**. Avoid multiple freeze-thaw cycles for XtenFect Reagent, Stock solution as it can lead to reduced transfection performances. If needed, the XtenFect Reagent, Stock solution can be aliquoted. The diluted XtenFect Reagent, Working solution could be **stored at 4°C** and should be used within **4 weeks** for optimal performance.

5. How to prepare ready-to-use XtenCHO™ Expression Medium?

The XtenCHO™ Expression Medium is **serum-free, animal-component-free, and chemically-defined**. The XtenCHO™ Expression Medium optimizes productivity with sustained balanced growth and viability post-transfection. The XtenCHO™ Expression Medium supports **small-and large-scale transient transfection** in XtenCHO™ Race Cells, enabling researchers to use the same medium for regular subculturing and for transient transfection.

The XtenCHO™ Expression Medium is formulated **without L-Glutamine and Anticlumping agent**. L-Glutamine is quickly hydrolyzed in medium, decreasing its actual concentration and generating byproducts of hydrolysis which can be toxic to the cells. To ensure optimal L-Glutamine concentration in medium, users should supplement the XtenCHO™ Expression Medium with **the recommended concentration of fresh L-Glutamine (8 mM)**.



The XtenCHO™ Expression Medium should also be supplemented with an **Anticlumping agent** (such as 0010057AE from Gibco) for subculturing prior to transfection and post-transfection. However, the Anticlumping agent **should not be added** in the XtenCHO™ Expression Medium **during the transfection** or it will result in decreased protein yields. It is highly recommended to remove the Anticlumping agent **during the last passage** before the transfection, by replacing complete XtenCHO™ Expression Medium (XtenCHO™ Expression Medium + 8 mM L-Glutamine + 0.5% Anticlumping agent) with XtenCHO™ Expression Medium **only supplemented with 8 mM L-Glutamine**.

6. How many passages are required before using XtenCHO™ Race Cells and how long can they be maintained in culture?

XtenCHO™ Race Cells should be passaged at least **three times** after thawing and should be cultured within the ranges specified in the XtenCHO™ Race Starter Kit User Guide, prior to transfection. Cells can be maintained up to **20 passages** without alterations in performance if maintained in accordance with the cell culture maintenance guidelines in the guide.

7. What is the expected doubling time for XtenCHO™ Race Cells?

The doubling time for XtenCHO™ Race Cells is about **21 hours**.

8. From which cell line are XtenCHO™ Race Cells derived?

XtenCHO™ Race Cell line is derived from a **CHO-K1 cell line**.

9. What are the normal growth characteristics of XtenCHO™ Race cells?

Within 2–3 passages after thawing, the doubling time of XtenCHO™ Race Cells should be of approximately **21 hours**. When cells are seeded at **0.2 – 0.3 x 10⁶ cells/mL**, viable cell density should be approximately **1 – 3 x 10⁶ cells/mL within 2 or 3 days**. Typically, cells should be subcultured each 2–3 days, when they reach a density of 2–3 x 10⁶ cells/mL. If cells are not growing within these approximate ranges, cell culture conditions will require further optimization. In such cases, different shaking speeds could be tested to determine which conditions allow optimal cell growth and high viability.



10. What is the expected viability for XtenCHO™ Race Cells and how can I analyse cell viability of XtenCHO™ Race Cells?

For both regular passaging and for transfection of XtenCHO™ Race Cells, **viability should be greater than 95%** to achieve best results with the XtenCHO™ Race Expression System. We recommend analyzing cell viability with trypan blue exclusion method.

11. For how many passages would you recommend using XtenCHO™ Race Cells before replacing them with new cells that were frozen at earlier passages?

Optimal and consistent expression results will be obtained when cells are cultured for less than **20 passages**.

12. What are the recommended culture conditions for XtenCHO™ Race Cells?

For cell line maintenance, we recommend passaging the cells every 2–3 days, at a cell density of **0.2 – 0.3 x 10⁶ cell/mL**. Cells should be passaged the day before transfection and seeded at a higher density (1.2 x 10⁶ cells/mL). On the day of transfection, cells are seeded at a density of **5 x 10⁶ cells/mL**. In case of cell clumping, let cell clumps settle in the vessel and only use the cells in suspension for subculturing. Indeed, the presence of clumps could alter significantly the accuracy of counting data. It is recommended to discard remaining cells as high-density cell cultures should not be used for seeding of maintenance flasks.

13. Can the XtenCHO™ Race Expression System be scaled down to smaller volumes?

Yes. Small scale transfections can be performed in vented **50 mL tubes** or in Deepwell plates with shaking. Reagents, plasmids and cells should be scaled down accordingly.



14. When subculturing XtenCHO™ Race Cells the day before a transfection, would you recommend using fresh medium to resuspend cells or just diluting cells to the desired cell density?

We recommend to change out the medium for XtenCHO™ Race Cells **the day before a transfection**. This is particularly important as the medium **should not contain Anticlumping agent** for transfection. On day of transfection, it is also advised to **use fresh medium**, without Anticlumping agent to resuspend cells. During the subculturing and amplification of the cells, it is also preferable to use fresh medium at each passage to reach optimal performance. However, to amplify large volume of cultures, **half of the medium volume** could be replaced for passages before the transfection, provided that Anticlumping agent is totally removed for the transfection.

15. What is the recommended density for the XtenCHO™ Race Cells at the time of transfection?

XtenCHO™ Race Cells should be subcultured at a density of **1.2 x 10⁶ cells/mL one day prior to transfection** to obtain a cell density of approximately **2–3 x 10⁶ cells/mL** on the day of transfection. These cells should be diluted to a final density of **5 x 10⁶ cells/mL** with fresh media and gently swirled to mix prior to transfection. It is recommended to **discard remaining cells** as high-density cell cultures should not be used for seeding of maintenance flasks.

16. Should I change the medium after the transfection?

No. The XtenCHO™ Race Expression System **does not require any medium exchange**. There is no need to remove transfection complexes or to change expression medium following transfection. No feed is required. However, the addition of the XtenCHO™ Race Enhancer, of XtenCHO™ Expression Medium two hours post-transfection, and the temperature **shift to 33°C on Day 1** after transfection, are necessary to achieve optimal results.



17. When should the XtenCHO™ Race Enhancer be added to the cultures? When should I add the XtenCHO™ Expression medium ?

The XtenCHO™ Expression Medium and the XtenCHO™ Race Enhancer should be added **two hours post-transfection** for best results. The XtenCHO™ Expression Medium and the XtenCHO™ Race Enhancer may be added to the flasks **with pre-warming at room temperature**. It is recommended to add first the XtenCHO™ Expression Medium and then the XtenCHO™ Race Enhancer in the transfected cultures. It is recommended to add the XtenCHO™ Race Enhancer closer to the timepoint for best performance.

18. Do I need to use two incubators at 37°C and 33°C or can I simply change the temperature of the 37°C incubator down to 33°C?

Optimal expression performances are obtained when using two incubators. However, for more convenience, the temperature can be lowered **from 37°C to 33°C** on the same incubator on **Day 1 post-transfection** (It can require some time to reach the desired temperature). The impact on the final expression yield will vary depending on the protein of interest.

19. Can I use a different medium for XtenCHO™ Race cells in my transfection?

The XtenCHO™ Race Starter Kit is optimized for use with the provided medium. This medium is **chemically-defined, animal component-free medium and designed to achieve high yields** of recombinant protein by transient transfection. Using another medium would not guarantee to achieve high rates of production for recombinant proteins, and could inhibit recombinant protein expression.

20. Can the amount of plasmid be changed to transfect XtenCHO™ Race Cells?

The protocol is optimized for use with the indicated DNA quantity **for optimal complexation with the XtenFect Reagent, Working solution**. However, it is possible to modify the quantity of DNA to transfect and conditions should be tested empirically for optimal results with your protein of interest.



21. For the XtenCHO™ Race Transfection protocol, what is the optimal time point for harvesting?

Harvesting time point depends on the nature of the protein expressed and the expression strategy (intracellular or secreted). Transfected cells can be kept until viability drops **under 50% or till Day 15 post transfection for flasks or Day 8 post transfection for plates**. Maintaining high cell viability at the time of harvest facilitates downstream processes, such as protein purification.

For stable secreted proteins, as recombinant antibodies, time of harvest occurs generally on **Day 13 to Day 15**. For intracellular or unstable proteins, it might be necessary to harvest culture at an earlier time point, to avoid toxicity or degradation. Optimal time to harvest should be determined empirically.

22. How should I process the cell supernatant?

It is recommended to harvest and centrifuge cultures for **5 min at 300 x g** to pellet cells. For recombinant protein whose expression is expected to be intracellular, **keep and freeze cell pellets at -80°C for further analysis**. Culture supernatants can be further clarified by centrifugation for **30 min at 5000-10000 x g** and/or filtration through a 0.22 µm filter. It is recommended to proceed quickly to purification of the protein of interest, or to freeze crude supernatants for later use.

23. How does the protein expression yield of the XtenCHO™ Race Expression System compare to ExpiCHO™ or CHO-S™ expression systems?

We observe **higher yields for the production of various recombinant proteins and therapeutic molecules** using the XtenCHO™ Race Expression System compared to CHO-S™ and ExpiCHO™ expression systems.



24. What amount of heavy and light chain encoding plasmids should be used to express recombinant antibody in XtenCHO™ Race Cells?

We observed that a **1:1 mass ratio of heavy- and light-chain encoding plasmids** can be generally used for expression of recombinant antibody with the XtenCHO™ Race Cells. Most of the cells will be **cotransfected** by the two plasmids.

However, the expression of the recombinant antibody may then be dependent on the expression rate of each chain and can be further optimized, if necessary, by modifying heavy to light chain plasmid ratio. We recommend cloning the heavy- and light-chain subunits separately **into the pXten1 expression plasmid** and before optimizing the ratios of the two plasmids.

25. Can XtenCHO™ Race Cells be used to generate Stable Cell lines?

Our XtenCHO™ Race Cells are derived from **CHO-K1 cells**, thus stable CHO-K1 selection methods are applicable to these cells. Note that XtenCHO™ Race Cells are engineered cells and are **resistant to geneticin**. It is recommended to use a different selection to establish stable cell lines.

26. Can a different transfection reagent be use for the transfection of XtenCHO™ Race Cells?

XtenFect Reagent, XtenCHO™ Race Enhancer and XtenCHO™ Expression Medium are optimized to work together to provide optimal protein expression levels. The XtenFect Reagent provides **high transfection efficiency of high-density cultures**. Combined with pXten1 Expression Vector, and transfected with the optimized protocol and reagents, **XtenCHO™ Race Cells sustain extended transient gene expression**, by enhancing plasmid-driven expression and allowing better plasmid maintenance.



27. Can the XtenCHO™ Expression Medium be used with other CHO Cell lines than the XtenCHO™ Race Cell line?

Yes, you may be able to adapt your CHO cells into XtenCHO™ Expression Medium. This change may increase productivity of your CHO cells and could allow **high-density growth**. However, there is no guarantee that CHO lines other than the XtenCHO™ Race Cell line will achieve the same levels of expression as the XtenCHO™ Race Cells. Indeed, XtenCHO™ Race Cells have been engineered to provide **optimal protein expression**. Combined with the pXten1 Expression Vector and transfected with the optimized protocol and reagents, **XtenCHO™ Race Cells sustain extended transient gene expression**, by enhancing plasmid-driven expression and allowing better plasmid maintenance.

28. How should I add reagents for the transfection?

High quality and purity plasmid DNA could be added all at once to the XtenCHO™ Race Cell culture. Then the prepared XtenFect Reagent, Working solution is added **drop-by-drop while swirling gently the flask**. Two hours post-transfection, the XtenCHO™ Expression **Medium, pre-warmed** to room temperature can be added all at once to the transfected cells. Similarly, the XtenCHO™ Race Enhancer pre-warmed to room temperature, should be added to the culture. It is recommended to add the XtenCHO™ Expression Medium **before** adding the XtenCHO™ Race Enhancer to the transfected cells. **Gently swirl** the flask to ensure even distribution.

29. Can I use minipreps to transfect XtenCHO™ Race Cells ?

We recommend using only **high quality, endotoxin-free plasmid** preparation for transfecting XtenCHO™ Race Cells as low quality and low purity plasmids will result in low transfection efficiency, lead to cell toxicity, and will dramatically affect bioproduction yield.



Ordering Information

Ordering information is provided below. The following reagents supplied in the XtenCHO™ Race Starter Kit and other products suitable for use with the kit are available separately. For more information, refer to proteogenix.science or contact the **Technical Support**.

Products	Amount	Reference
XtenCHO™ Race Starter Kit, includes: <ul style="list-style-type: none">■ XtenCHO™ Race Cells■ XtenCHO™ Expression Medium■ XtenFect Reagent, Stock solution (to dilute prior to use)■ XtenFect Reagent, Working solution (to reconstitute prior to use)■ XtenCHO™ Race Enhancer■ pXten1 Expression Vector■ Antibody Expression Positive Control Vector	2 x 1 mL 1 x 1 L 1 x 0.51 mL 3 x 1.44 mL 3 x 1.2 mL 1 x 10 µL 1 x 150 µL	PX-XTE-R-001
XtenCHO™ Expression Medium	1 x 1 L	PX-XTE-002
XtenCHO™ Race Transfection Kit, includes: <ul style="list-style-type: none">■ XtenFect Reagent, Stock solution (to dilute prior to use)■ XtenFect Reagent, Working solution (to reconstitute prior to use)■ XtenCHO™ Race Enhancer	1 x 0.51 mL 3 x 1.44 mL 3 x 1.2 mL	PX-XTE-R-003
XtenCHO™ Race Cells	1 x 1 mL	PX-XTE-R-004
XtenCHO™ Race Cell Bank Pack	6 x 1 mL	PX-XTE-R-005
XtenCHO™ Expression Medium Pack	10 x 1 L	PX-XTE-006



+33 (0)3 90 20 54 70

support-xtencho@proteogenix.fr

www.proteogenix.science

